# BLUE LIFT

Polyp Lifting Solution 10 ml

REF BL1

STERILE

## 1 - NAME AND ADDRESS OF THE DISTRIBUTOR

DERM TECH FRANCE

210 Avenue Aristide Briand  
93320 Pavillons-sous-bois - France  
Tel : +33 1 48 02 61 63 – Fax : +33 1 48 02 61 71  
info@derm-tech.com

## 2 - IDENTIFICATION OF THE PRODUCT

### a-Name

Blue Lift - Polyp Lifting Solution 10 mL sterile medical device

### b-Composition

Sodium Hyaluronate, Indigo carmine, Saline.

### c-Package Contents

REF BL1 : Solution in a clear glass vial for antibiotics containing 10 ml (filled to 9.7ml), sealed with a butyl rubber stopper and closed with an aluminum cap (crimp vial). Blue Lift is sterilized and filled into a ready to use and disposable vial Packaged in a box containing 10 single use vials.

### d- Name and address of the manufacturer

#### DERM TECH FRANCE

210 Avenue Aristide Briand – 93320 Pavillons-sous-bois - France  
Tel : +33 1 48 02 61 63 – Fax : +33 1 48 02 61 71

info@derm-tech.com

### e- Device classification

Class IIa medical device sterile

### f- Method of sterilization

Sterile product for single use only. Steam Sterilization Cycle prion, 134 °C for 18 minutes.

## 3 - INTENDED USE

### When to use this device (indications for use)

Blue Lift is a transient use product which is injected into the submucosa beneath the lesion to be excised. Blue Lift is effective immediately after injection.

Following injection into the submucosa Blue Lift pushes the mucosa away from the muscularis propria (muscle) layer and enables a snare to be positioned around the base of the polyp.

The cushion of Blue Lift within the mucosa protects the underlying muscle from thermal damage, which can be caused by heat generated during the electro cauterization of the polyp.

The viscous Blue Lift solution also acts as a tamponade, constricting blood vessels and thereby minimizing bleeding from the wound. It is recommended that any marginal tissue surrounding the polyp or lesion extraction site is removed.

## 4- WARNING!

### a- When not to use this device: known contraindications

#### This device SHOULD NOT BE USED:

- Patients with known sensitivity to hyaluronate or indigo carmine.
- By an individual not trained in preoperative marking techniques of lesions or tumors or radiolabeling

### b -Precautions for use

Before using the device, check the physical appearance of the suspension (packaging has no visible defect, colour and texture of liquid and homogeneity after shaking).

The Endoscopist injecting the Blue Lift must be experienced in the technique of submucosal injection for resection of superficial gastrointestinal lesions.

Sodium hyaluronate is manufactured by fermentation of Streptococcus equi and rigorously purified. However, the physician should consider the immunological and other potential risks that can be associated with the injection of any biological material.

The safety of Blue Lift has not been established in pregnant or lactating women or children under 18 years of age.

Needles should be sterile and are not reused because Blue Lift may become turbid due to precipitation in the presence of cationic agents such as benzalkonium chloride or detergents residual in the needle following re-sterilization.

Blue Lift must not be injected into blood vessels because it has the potential to occlude the vessels, which could result in an embolism or infarction.

### c- Warnings and Precautions

Product quality and safety cannot be guaranteed if the product has been altered.

This solution is for single use only. Any open product must be disposed of after use.

**IF YOU HAVE ANY FURTHER QUESTIONS, CONTACT YOUR DOCTOR OR PHARMACIST FOR ADVICE.**

## 5- HOW TO USE THIS DEVICE

### a-Instructions for use

Blue Lift is a solution of sodium hyaluronate with Indigo Carmine in a saline solution contained in a pre-filled vial for submucosal injection to expand, deepen and highlight the submucosal layer beneath the tissue to be excised.

Blue Lift can be injected via a 23GA Sterile Endoscopic Needle and through an endoscope. Blue Lift is supplied in a 10ml vial which can be drawn into a syringe.

The sodium hyaluronate contained in Blue Lift is manufactured by biosynthesis using the natural bacterium *Streptococcus equi*. The polysaccharide is not modified in any way. Physicochemical tests have confirmed that the molecular composition and configuration is the same as that of sodium hyaluronate found naturally in the body

### b-Dosage and Administration:

Blue Lift is injected via an endoscope using a sterile endoscopic needle (endoscopic needle injector 21GA, 23GA or 25GA)

Position the endoscope, with a biopsy channel of sufficient diameter to accommodate the injection catheter and needle, close to the lesion to be resected in the patient.

Draw the Blue Lift into a suitably sized syringe.

Hold the syringe upright and allow all air bubbles to rise to the end of the syringe.

Gently expel the air from the syringe.

Securely attach the syringe to the injection needle. Prime the injection needle catheter with Blue Lift and expel any air.

Introduce the sterile injection needle through the endoscope.

Inject Blue Lift into the submucosal layer beneath the tissue to be excised in one or more injections.

Consumables have been in contact with the patient must be discarded with DASRI\* in compliance with the legislation.

(\*) DASRI: Care Activity Waste for Infectious risk.

Blue Lift bottle can be disposed of as normal waste.

## 6- UNDESIRABLE EFFECTS

### Adverse reactions:

Rarely, an inflammatory reaction could occur which may or may not be caused or associated with Blue Lift.

**REPORT ANY UNDESIRABLE EFFECT NOT MENTIONED IN THIS LEAFLET TO DERM TECH FRANCE**

## 7- STORAGE AND SHELF LIFE

### a-Shelf life

DO NOT USE AFTER THE EXPIRY DATE PRINTED ON THE PACKAGING. SOLUTION FOR SINGLE USE ONLY. ANY OPENED PRODUCT MUST BE DISPOSED OF AFTER USE.

### b- Special precautions for storage

Store at a temperature between 0°C and 4°C in a cool place, dry place protected from Direct Sunlight. Do not Freeze.

This is a light-sensitive product which must be stored protected from light. Please adhere to storage instructions. It is possible that incorrect storage could impact on product coloration.

### c- Single use product

This product is for single use only. Do not reuse, re-sterilize or repackage

### d-Storage after opening

This product is for single use only. Any opened product must be discarded after use and may not be stored for later use. Re-use of a product after opening may cause a risk of infection

### e- Warning in case of visible signs of deterioration

In case of deterioration of the vial, product quality and safety cannot be ensured. Before using, check the integrity of the vial


## 8- GENERAL CONDITIONS

This device is for single use only. Any use not indicated in the leaflet is not guaranteed by the manufacturer.

## 9- LEGEND OF SYMBOLS USED

The symbols are those used in the EEC EN ISO 15223-1 standard

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