
CENDER'S - Endoscopic marker 5mL

REF

JF/185 CENDER

1- NAME AND ADDRESS OF THE MANUFACTURER

DERM TECH FRANCE

DERM TECH FRANCE

210 avenue Aristide Briand, 93320 Les Pavillons-sous-Bois - France
Tel.: +33 1 48 02 61 63 - Fax: +33 1 48 02 61 71
info@derm-tech.com

2- IDENTIFICATION OF THE PRODUCT

a- Name

Cender's endoscopic marker 5mL – Sterile pigment for endopigmentation

b- Composition

Water for Injections, Peat charcoal (in compliance with the US Food Chemicals codex).

c- Package Contents

10% peat charcoal suspension in Water for Injections.

REF.: JF/185CENDER: 10% suspension in a clear glass vial for antibiotics containing 8ml (filled to 5ml), sealed with a butyl rubber stopper and closed with an aluminium cap (crimp vial).

Packaged in a box containing 10 vials.

d- Device classification

Class IIb medical device (in accordance with rule 8 of Annex IX of EEC Directive 93/42).

e- Method of sterilisation

Sterile product. Steam sterilisation, prion cycle, 134 °C for 18 minutes.

3- INTENDED USE

When to use this device (indications for use)

This medical device is a pigment for endopigmentation and has the advantage of not diffusing into the surrounding tissues.

a- Breast tumour marking prior to surgery and chemotherapy.

This technique allows the pathologist, during frozen section examinations, to guide the surgeon when there is no longer any residual tumour.

b- In gastroenterology, endoscopic tattooing of polyps or tumours for surgical identification and/or long-term endoscopic and histological follow-up after the excision of polyps or tumours.

c- Thyroid marking

4- WARNING!

a- When not to use this device: known contraindications

This device SHOULD NOT BE USED:

- if there is a history of allergy to one of the components
- by an individual not trained in preoperative marking techniques of lesions or tumours or radiolabeling.

b - Precautions for use

Before using the device, check the physical appearance of the suspension (packaging has no visible defect, colour and texture of liquid and homogeneity after shaking).

c- Warnings and Precautions

Sterility cannot be guaranteed if the contents have been altered.

This solution is for single use only. Any vial or syringe that has been opened must be disposed of as the product is no longer sterile.

IF YOU HAVE ANY QUESTIONS, CONTACT YOUR DOCTOR OR PHARMACIST FOR ADVICE.

d- The product comes in a sterile condition. The product is a single use.

It must in no case be reused to avoid patient contamination.

(Cross contamination or loss of sterility)

5- HOW TO USE THIS DEVICE

a- Accessories to be used with this device

The essential accessories are: the needle (it should be delivered sterile in double packaging). For each accessory, See the individual instructions for use for each accessory.

b-Instructions for use:

REF: JF/185CENDER:

- Shake the vial for 1 minute before using.
- Remove the vial from its double packaging and remove the cap from the vial until the centre of the stopper is visible. Be careful to keep your hand clear so as not to cut yourself.
- Using a sterile syringe inserted in the centre of the stopper, withdraw the required quantity of suspension.

b.1- For marking breast tumours prior to surgery and chemotherapy, 1 ml of the suspension is injected into the tumour or around the tumour nodule during the initial biopsy for diagnostic purposes. This is required so that the tumour can be observed following breast cancer treatment by neoadjuvant chemotherapy.

b.2- For colonic tattooing (polyps or tumours)

The lesion perimeter is marked at 4 points by the submucosal endoscopic injection of 0.5 to 1 ml of 10% suspension. The injection of the suspension makes it possible to have long-term endoscopic and histological follow-up by colonoscopy following tumourectomy or polypectomy and to guide the surgeon by tattooing the lesion.

b.3- Use the product without adding any other substance

b.4- What to do if the product becomes blocked:

- Remove the needle from the mucosa (the catheter can remain in the endoscope).
- Aspirate a small amount of air into the syringe.
- Shake the syringe firmly.
- Expel all the air from syringe and catheter until a drop of product appears at the needle tip.
- Resume the marking.

Consumables have been in contact with the patient must be discarded with DASRI* in compliance with the legislation.

(*) DASRI: Care Activity Waste for Infectious risk.

Cender's bottle can be disposed of as normal waste.

6- UNDESIRABLE EFFECTS

Moderate risk of pain during the intratumoural injection.

REPORT ANY UNDESIRABLE EFFECT NOT MENTIONED IN THIS LEAFLET TO DERM TECH FRANCE.

7- STORAGE AND SHELF LIFE

a- Shelf life

DO NOT USE AFTER THE EXPIRY DATE PRINTED ON THE PACKAGING.

b- Special precautions for storage

Store at a temperature between 5°C and 45°C in a clean, dry place protected from light.

c- Single use product

Single use product.

Do not reuse, do not resterilise, do not repackage.

d- Storage after opening

Solution for single use. Any opened vial or syringe must be discarded after use and may not be stored for later use.

e- Warning in case of visible signs of deterioration

In case of deterioration of the vial or syringe, sterility cannot be ensured. Before using, check the integrity of the vial or syringe.

8- GENERAL CONDITIONS

This device is for single use only.

Any use not indicated in the leaflet is not guaranteed by the manufacturer.

9- LEGEND OF SYMBOLS USED

The symbols are those used in the ISO15223-1.

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