



TECHNICAL DATA SHEET

BLUE LIFT

Endoscopic Marker



➤ MANUFACTURER

DERM TECH France

210 Avenue Aristide Briand, 93320 Les Pavillons-sous-Bois

Tel: +33 (0)1 48 02 61 63 / Fax +33 (0)1 48 02 61 71

E-mail: nohaman@derm-tech.comaicha@derm-tech.com

Medical Device Vigilance Coordinator

General Manager

Tel: +33 (0)1 48 02 61 63

Fax: +33 (0)1 48 02 61 71

E-mail: nohaman@derm-tech.com

➤ INFORMATION ABOUT THE PRODUCT

Trade name**DEEP BLUE**

Endoscopic marker – 10 mL

Polyp Lifting Solution

Sodium Hyaluronate, Indigo carmine, Saline.

Indications :

Deep Blue is a transient use product which is injected into the submucosa beneath the lesion to be excised. Deep Blue is effective immediately after injection.

Following injection into the submucosa Deep Blue pushes the mucosa away from the muscularis propria (muscle) layer and enables a snare to be positioned around the base of the polyp.

The cushion of Deep Blue within the mucosa protects the underlying muscle from thermal damage, which can be caused by heat generated during the electro cauterisation of the polyp.

The viscous Deep Blue solution also acts as a tamponade, constricting blood vessels and thereby minimising bleeding from the wound. It is recommended that any marginal tissue surrounding the polyp or lesion extraction site is removed.

Reference:

REF BL

Product description

10 BLUE LIFT VIALS

Packaging sizes**Individual packaging:**

Antibiotic glass vial, 10 mL, in a box

Packaging sizes:

Cardboard box containing 10 Vials

Order unit: Only boxes of 10 units**Storage:** Controlled temperature: (0°C - 4°C) in a cool place

TECHNICAL DATA SHEET

BLUE LIFT

Protect from direct sunlight and humidity

Shelf life: 12 months

Class I sterile medical device – Manufactured by Derm Tech France (France) – CE marking granted by Eurofins, Notified Body no.: 0477 – Safety in use, recommendations for use (instructions for use, precautions and contraindications): refer to the information leaflet – Device not reimbursed by healthcare insurance policies.

1/2

Composition:

Sodium Hyaluronate, Indigo carmine, Saline

Presence of:

Latex:	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No
Product of animal origin:	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No
PVC:	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No
Phthalate and DEHP:	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No
Bisphenol A:	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No

Sterility:

Sterile: Yes No



Sterilisation method: Autoclave
 Number of uses: Single Use



Storage:

Storage conditions: Controlled temperature : Store at a temperature between 0°C and 4°C in a cool place, dry place protected from Direct Sunlight. Do not Freeze.

Protect from direct sunlight and humidity

Shelf life: 12 months

Other data:

GMDN Classification:	available	
Bar code present:	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No
Reference standard used:	<input type="checkbox"/> HIBC	<input checked="" type="checkbox"/> GS1-128

With information:

Reference:	<input checked="" type="checkbox"/> Yes	<input type="checkbox"/> No
Batch no.:	<input checked="" type="checkbox"/> Yes	<input type="checkbox"/> No
Expiry date:	<input checked="" type="checkbox"/> Yes	<input type="checkbox"/> No

Use:



Safety in use: Refer to Information Leaflet

Recommendations for use (instructions for use, precautions and contraindications): Refer to Information Leaflet

Confidential document

Application date: 17/03/2020

Technical data sheet - Derm Tech – DEEP BLUE



TECHNICAL DATA SHEET

BLUE LIFT

Indications (CE marking intended use): Refer to Information Leaflet

Manufacturer:

Manufacturer: DERM TECH France

Country: France

CE marking

CE marking granted by Eurofins, Notified Body no.:0477 According to European directive 93/42/EEC and its amendments. Class IIa sterile medical device.