


TECHNICAL SHEET

Endoscopic marker – Cender's



➤ MANUFACTURER

DERM TECH France 210 avenue Aristide Briand – 93320 Les Pavillons-sous-Bois Tél : 01 48 02 61 63 / Fax 01 48 02 61 71 E-mails : nohaman@derm-tech.com aicha@derm-tech.com	Quality Representative Quality manager Tél : 01 48 02 61 63 Fax : 01 48 02 61 71 E-mails : nohaman@derm-tech.com
--	---

➤ PRODUCT INFORMATION

Commercial name :

CENDER'S
Endoscopic Tattoo – 5 ml
For injection for endopigmentation
10 % charcoal suspension in sterilized water

Indications: This medical device is a pigment for endopigmentation, the interest of which is not to diffuse in the surrounding tissues.

- In gastroenterology endoscopy labeling of polyps and tumors for a surgical localization with or without a prolonged endoscopy and histology follow-up after polyps or tumor exeresis.

Presurgery and prechemiotherapy labeling of breast tumors. This technique allows the anatomopathologist, when doing an extemporaneous examen, to guide the extend of the surgery when there is no more residual tumoral nodule.

Référence :	Description du produit	Conditionnement
JF/185 CENDER'S	10 Bottles of cender's	Individual packaging: Antibiotic bottle 5 ML Overall packaging box of 10 bottles

Minimum Delivered Quantity : 10

1/2

TECHNICAL SHEET

Composition :

Micronized charcoal suspension , PPI water (fit with US Food chemicals codex)

Is there any :

Latex :	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No
Animal Origin or Residues :	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No
PVC :	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No
Phtalate and DEHP :	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No
Bisphénol A :	<input type="checkbox"/> Yes	<input checked="" type="checkbox"/> No



Sterility :

Steril : Yes No



Nombre d'utilisations : Single use

Storage conditions :

Conditions of conservation: Store in a controlled room temperature (5-45°C)
In a dry and safe light place

Validity : 3 years

Other informations :

GMDN code : available

Barcode data matrix : Yes No



Using Advices: Cf. Instruction manual — Just before using, shake the product for more than 1 minute in order to mix the micronized particles. Inject Sterimark through a 0.7 mm injection needle

Indications according to CE mark): Cf Instruction manual — digestive Tumours, polyps or mucosa tattooing, for endoscopy, surgery or histology follow up.

Warnings : Cf. Instruction manual — This product must not be in case of previous allergy to one of the components, and by a person not technically formed to preoperatory labelling techniques of lesions or tumours, or to radio labelling.

Manufacturer: Derm Tech France
Manufacturer: France
Country:



TECHNICAL SHEET

Cender's

CE Mark:

CE Mark, by Notified body EUROFINS,
Number: 0477

Certificate No EPT 0477 MDD 21/4238